

dataLink pro

Petri dish labeling solution

Réf. 439 050



In compliance with:

21 CFR
Part 11

Traceability:



Manufactured:



The dataLink pro pack facilitates Petri dish identification for optimum traceability at all steps of microbiological analysis. The dataLink software controls and transmits data on a label. Semi-automated labelling on the edge of the dish is great for large production analyses.

OVERVIEW

- Sample identification and traceability
- Automation of counting settings
- Import of raw data
- Customized labeling
- Semi-automatic labeling tray

SPECIFICATIONS

- Type of labeling: Semi-Automated
- Number of labels applied per minute: 20
- Connection: Network or USB
- 1 Hi Speed USB 2.0 port, 1 front USB port
- 1-year warranty
- In compliance with 21 CFR Part 11
- GLP-compliant use (Good Laboratory Practice)
- Manufacturing : CE, UKCA

SOFTWARE

- Data types: Volume, dilution, hour, plating mode, user
- Data formats:
Input: alphanumerical
Output: alphanumerical, barcode, datamatrix

PARIS

Phone: +33 1 34 62 62 61
info@interscience.com

FRANKFURT

Phone: +49 611 7238 7770
sales.germany@interscience.com

BOSTON

Phone: +1 781 937 0007
sales.usa@intersciencelab.com

SHANGHAI

电话: +86 178 2123 6642
sales.china@interscience.cn

SINGAPORE

Phone: +65 6977 7232
sales.asia@interscience.com

TOKYO

Phone: +81 3 6712 9715
sales.japan@interscience.com

- Data Import/Export:
Automatic: csv
Manual: txt, xlsx
- Language: English
- User account management: Local or LDAP
- Audit trail

PRINTER

- Printing type: Thermal transfer
- LCD screen
- Label dimensions: 50 x 5 mm
- Number of characters: Maximum of 73 alphanumeric characters
(varies depending on special characters/uppercase letters/barcode types)
- Printing resolution: 600 dpi
- Voltage – Frequency: 100-240 V~ 50-60 Hz
- Max Power : 200 W

PC MINIMUM REQUIREMENTS

- Operating system: Windows™ 10 or 11 (or higher)
- Processor: Intel i5, 2,8 GHz, quad Core (or higher)
- RAM: 4 GB
- Equipment: 4 free USB ports (USB hub compatible) and/or 1 RJ45 port
- Screen resolution: 1280 x 1024 pixels (or higher)

WEIGHT AND DIMENSIONS

- Dimensions (w x d x h): 30 x 56.9 x 25.5 cm, net weight: 12 kg
- Box (w x d x h): 70 x 35 x 40 cm, gross weight: 16 kg

PACK CONTENT

- 1 dataLink pro printer
- dataLink software
- 1 power cord
- 1 USB cable
- 2 user manuals
- 1 thermal transfer ribbon
- 1 chuck
- 3 rolls of 2 500 adhesive labels

OPTIONS / ACCESSORIES

- Ref. 439 051 **Stickers for printer (dataLink pro):** 3 rolls of 2 500 labels
for dataLink pro printer - Label size: 50 x 5 mm - Quantity: 3 x 2500
- Ref. 439 052 **Film for printer (dataLink pro):** Thermal transfer film for dataLink pro printer
- Ref. 439 170 **Petri dish automatic identification kit:**
Bar-code (1D) and datamatrix (2D) reader with a turntable for automatic identification of Petri dishes
with ScanStation
- Ref. 522 000 **Barcode reader:** Barcode reader (1D / 2D)

PARIS

Phone: +33 1 34 62 62 61
info@interscience.com

FRANKFURT

Phone: +49 611 7238 7770
sales.germany@interscience.com

BOSTON

Phone: +1 781 937 0007
sales.usa@intersciencelab.com

SHANGHAI

电话: +86 178 2123 6642
sales.china@interscience.cn

SINGAPORE

Phone: +65 6977 7232
sales.asia@interscience.com

TOKYO

Phone: +81 3 6712 9715
sales.japan@interscience.com