

BagTips

Straw tips for pipet

Ref. 252 020 (BagTips Protect) - Ref. 252 019 (BagTips Regular) - Ref. 252 024 (BagTips Jumbo)



BagTips straw tips, used with a pipet, make it easier to collect samples by effectively reaching the bottom of the bag, while ensuring optimal safety. Available in three versions (Protect, Regular, and Jumbo), they are suitable for all types of liquid samples and containers. They ensure the reliability of the analyses. BagTips Protect, with an integrated filter, effectively protects the pipet from aerosol effects and cross-contamination.

SPECIFICATIONS

- Compatible with most straw pipets
- Container types: Suitable for sampling from any type of container
- Sample types: Liquid samples
- Material: Polypropylene
- Gamma-ray treatment: 5 to 10 kGy
- Storage conditions: In a cool, dry place, away from moisture and light
- Complying with ISO 7218, ISO 6887-1, ISO 4833-2 et FDA BAM (Bacteriological Analytical Manual)
- Compliant-use in accordance with Good Laboratory Practice (GLP)

BagTips Protect

Ref.: 252 020

- Integrated filter
- Nominal volume: 2 mL
- Volume range⁽¹⁾: 100 to 2 000 µL
- Filter porosity: 40 µm
- Hydrophobic filter
- Length: 190 mm ± 3 mm
- Box (w x d x h): 27.5 x 22.5 x 11 cm
- Gross weight: 1.24 kg
- Quantity per pouch: 25
- Quantity per box: 1 000

(1) varies depending on the pipet

FT_005686_Technical File BagTips_EN_01-02 Information and pictures are not binding. INTERSCIENCE reserves the right to change or improve his specifications without any notice. Please consult the website www.interscience.com for more information.
RCS 950 356 220 Versailles. INTERSCIENCE SAS F78860

BagTips Regular

Ref.: 252 019

- Nominal volume: 2 mL
- Volume range⁽¹⁾: 100 to 2 000 µL
- Length: 190 mm ± 3 mm
- Box (w x d x h): 27.5 x 22.5 x 11 cm
- Gross weight: 1.24 kg
- Quantity per pouch: 25
- Quantity per box: 1 000

BagTips Jumbo

Ref.: 252 024

- Nominal volume: 2.5 mL
- Volume range⁽¹⁾: 100 to 2 500 µL
- Length: 240 mm ± 3 mm
- Box (w x d x h): 27.5 x 22.5 x 11 cm
- Gross weight: 1.20 kg
- Quantity per pouch: 25
- Quantity per box: 1 000

APPLICATIONS

	BagTips Regular	BagTips Jumbo	BagTips Protect
Barrier against droplets, aerosols, and liquids	-	-	✓
Elimination of the cross-contamination risk from the pipet	-	-	✓
Quality indicator	● ● ●	● ● ●	● ● ●
Pathogen Analyses	● ● ●	● ● ●	● ● ●
400 mL-Bag	● ● ●	● ● ●	● ● ●
2 000/3 500 mL-Bag	● ● ●	● ● ●	● ● ●

ASSOCIATED PRODUCTS

Reference	Title	Description
251 091	BagPipet	Pipet with fixed volume

(1) varies depending on the pipet

FT_005686_Technical File BagTips_EN_01-02 Information and pictures are not binding. INTERSCIENCE reserves the right to change or improve his specifications without any notice. Please consult the website www.interscience.com for more information.

RCS 950 356 220 Versailles. INTERSCIENCE SAS F78860